

## Experiment 2 – Graphing Density

### Discussion

Density measures the “compactness” of material. For example, lead has a high density, and Styrofoam has a low density. Mathematically, this compactness is expressed as mass per unit volume. In chemistry, we use grams and milliliters:  $d = \text{g/mL}$ . Density is an intrinsic value; it does not depend on the amount of sample taken. We will take advantage of this by measuring the density of various sample sizes and averaging their densities.

### Procedure

1. Obtain 5 pieces of the same object and record its name on the data sheet.
2. Weigh each object and record the values on the data sheet.
3. Choose the appropriate size graduated cylinder (smallest size that will hold the object plus enough water). Add enough water to the graduated cylinder be able to cover your largest sample. Record the volume to the highest precision (0.1 mL or better).
4. Carefully add the sample to the graduated cylinder. There are two things to watch out for: breaking the cylinder and splashing water out. Tilting the cylinder and gently sliding the object in minimize both of these risks. Record the new volume.
5. Repeat with each sample piece. If the sample pieces are small, the pieces can remain in the graduated cylinder until all sample pieces have been added to the cylinder.
6. Determine the density of each piece and the average density.
7. Using Microsoft Excel®, graph the cumulative data. Use the largest values of mass and volume to determine your x and y scales.

### Excel® Procedure

Note that various versions of Excel® may function a bit differently from the directions outlined below (which work on department-owned laptop computers):




Figure 1. Creating a Scatter plot in Excel®.

Name: \_\_\_\_\_

Section: \_\_\_\_\_

Put the title for your x-axis (include units) in one Excel® cell (box). In the cell to the right, put the title for your y-axis. Using these boxes as headings, input the numeric data (like a table) in the cells under these titles (each box should contain one number; each row represents one data point in x,y format). Click and drag your mouse to highlight just the numeric boxes. From the “Insert” tab, choose a “Scatter” plot. (See example, below.)

Your graph must include a meaningful Chart Title and Axis Titles (with units). These Chart Elements can be added to your graph by clicking on the “+” icon in the upper right corner of your graph. Your instructor may request additional Chart Elements.

To add a Trendline, right click on any data point on your graph and choose “Display Trendline” from the menu that appears. The format trendline pane will appear on the right side of your screen. Linear should be selected by default. From this pane, you should check the box next to “Display Equation on chart.” Your instructor may also ask you to check the box for “Display R-squared value on chart.”

**Data and Calculations for Experiment 2**

Data Sheet for Density of an Element sample

Name of Element: \_\_\_\_\_

*Using Microsoft Excel<sup>®</sup>,  
graph the following:*

| Sample # | Object Mass (g) | Initial mL H <sub>2</sub> O | mL H <sub>2</sub> O w/ Object | Volume object (mL) | Density (g / mL) | Cumulative Sample #s | Cumulative volume (mL) (x-axis) | Cumulative object mass (g) (y-axis) |
|----------|-----------------|-----------------------------|-------------------------------|--------------------|------------------|----------------------|---------------------------------|-------------------------------------|
| 1        | _____           | _____                       | _____                         | _____              | _____            | 1                    | _____                           | _____                               |
| 2        | _____           | _____                       | _____                         | _____              | _____            | 1 + 2                | _____                           | _____                               |
| 3        | _____           | _____                       | _____                         | _____              | _____            | 1 + 2 + 3            | _____                           | _____                               |
| 4        | _____           | _____                       | _____                         | _____              | _____            | 1 + 2 + 3 + 4        | _____                           | _____                               |
| 5        | _____           | _____                       | _____                         | _____              | _____            | 1 + 2 + 3 + 4 + 5    | _____                           | _____                               |

Average Density from Table = \_\_\_\_\_

Average Density from Graph (slope of line) = \_\_\_\_\_

Be sure to show your properly formatted graph on the laptop or your device to your instructor to receive credit for this part of the experiment OR they may instruct you to take a photo of your graph to submit with this report to Canvas.

**Questions**

1. Which would work better in this experiment as an unknown solid whose density is to be determined, wood chips or small quartz rocks? Explain your choice.
2. Why is it best to use a smaller graduated cylinder as opposed to a larger graduated cylinder for this experiment?
3. How well does the average density from the table and density from the slope of the graph compare? Which value is closer to the accepted density of your metal in the table below (from the *CRC Handbook of Chemistry and Physics*)? Calculate the percent error between your better value and the accepted value.

| Element      | Accepted Density (g/mL) |
|--------------|-------------------------|
| Aluminum, Al | 2.70                    |
| Copper, Cu   | 8.96                    |
| Lead, Pb     | 11.34                   |
| Zinc, Zn     | 7.14                    |

4. What is the density of a 9.343 gram piece of metal that causes the level of water in a graduated cylinder to rise from 5.1 to 8.1 mL when the metal is submerged in the water? Consider significant figures when doing the calculation.